



TFS Mondo Re Circulating Chemical Portable Toilet



This high quality Re circulating Chemical Toilet has been a hit for years in Europe and is now taking the rest of the world by storm. Made from best materials only, it complies with the highest hygiene and comfort requirements and is built to last with low maintenance and total reliability. This Mondo Re Circulating Toilet includes:

Specifications:

Toilet seat • Lid, • Waste Compartment • Anchor points for secure lifting • larger waste capacity to hold more waste, therefore fewer cleans needed (241 Litres) • Re Circulating Hand Pump • Extra Strengthened door hinges • Double toilet paper holder • Please note – Urinals are not favoured in the UK market, these have no fresh water to flush urine through, however these are available on request at an additional charge • Please note if you opt for a urinal a sink cannot be fitted.

HOW TO WORK A CHEMICAL TOILET: Once you have received your chemical toilet, there are a few simple start up instructions to prepare for use: This unit works on a re circulating basis, meaning the fluid in the base of the waste tank is recycled and used for flushing, to set up your toilet you will need to: Lift the toilet seat lid completely and expose the base of the tank, you will see at the bottom of the tank a hose and filter, which is attached to the pull flush handle. Take your chemical additive and pour into the base of the tank the required amount, if using chemical supplied by us, you will see a small section on the side of the bottle, squeeze into that section the required amount (usually between 25-35ml) then pour into the base of the tank, if using -an alternative chemical, follow the guidelines on that – once the chemical has been added, you will need 8-10 litres of water (approx.) to cover the filter bag attached to the hose/handle situated in the base of the tank. On completion of the above, place the tank lid down and pull the hand pump flush handle to expel any air and wait for the water/chemical to pull up through the unit to flush, you can continue to use this until the unit becomes full, regular checks on the tank can be made by exposing the tank, you will notice the handle will become reluctant to flush if the unit becomes too full.

FOR EMPTYING:

For the disposal of waste, we recommend you contact a “Wet Waste Disposal” company within your area – these companies can be found in a local directory or via a Google search.

Unit Dimensions (HxWxD) 2.33m x 1.22m x 1.20m

Weight: 78 kgs

OPTIONAL EXTRAS

Bulkhead Dispenser (Can be filled with any hand soap or hand sanitizer gel)

Wall Mounted Deb Sanitizer Dispenser with 1 Litre Foam Refill

1 Litre of Chemical

Cold Water Foot Pumped Sink

Hot Water Foot Pumped Sink